

MODULAR EXTERNAL FIXATOR SYSTEM

For Common Use

Product Features

- With the characteristic of advanced technology and convenient usage.
- Make the frame compactly.
- Simple and effective usage.
- Less tools, effectively reduce the space occupancy.
- Simple structure, so that doctors can get more freedom.



MODEL- EXFIX XV

Modular External Fixator System (For Common Use)-Type III


REF	SIZE	IMAGE
MOI 55031084	Pin to Rod Coupling $\Phi 8 \times 4$ mm	
MOI 55038115	Pin to Rod Coupling $\Phi 11 \times 5$ mm	
MOI 55032088	Ringent Rod to Rod Coupling $\Phi 8 \times 8$ mm	
MOI 55039011	Ringent Rod to Rod Coupling $\Phi 11 \times 11$ mm	
MOI 55033845	Ringent Pin to Rod Coupling $\Phi 8 \times 4-5$ mm	
MOI 55040056	Ringent Pin to Rod Coupling $\Phi 11 \times 5-6$ mm	
MOI 55035088	Rod to Rod Coupling $\Phi 8 \times 8$ mm	
MOI 55042011	Rod to Rod Coupling $\Phi 11 \times 11$ mm	
MOI 55034084	T-clamp $\Phi 8 \times 4$ mm	
MOI 55041115	T-clamp $\Phi 11 \times 5$ mm	

MODEL- EXFIX XV


Modular External Fixator System (For Common Use)-Type III

REF	SIZE	IMAGE
MOI 15010011	Universal Clamp Closed	
MOI 15013011	Universal Clamp Open	
MOI 15011011	UNIVERSAL CLAMP	
MOI 15012011	TRANSVERSAL CLAMP	
MOI 15009011	TUBE TO TUBE CLAMP	
MOI 55321008	Semi-circular Trapezoidal Curved Rod Φ 8 mm	
MOI 55321011	Semi-circular Trapezoidal Curved Rod Φ 11 mm	
MOI 55322008	Trapezoidal Curved Rod Φ 8 mm	
MOI 55322011	Trapezoidal Curved Rod Φ 11 mm	
MOI 55323008	Curved Rod Φ 8 mm	
MOI 55323011	Curved Rod Φ 11 mm	

Recommended Frame Configurations




DESCRIPTION	REF	QTY
Standard Bone Screws	MOI 15417130	2
T-clamp	MOI 55034084	1
Standard Bone Screws	MOI 15418130	3
Pin to Rod Coupling	MOI 55031084	3
Connecting Rod	MOI 65340300	1
Used for: Tibia Comminuted		



DESCRIPTION	REF	QTY
Standard Bone Screws	MOI 15417130	4
Pin to Rod Coupling	MOI 55038115	4
Connecting Rod	MOI 65340150	1
Rod to Rod Coupling	MOI 55042011	2
Connecting Rod	MOI 65340200	2
Used for : Femur, Tibia		

Recommended Frame Configurations



DESCRIPTION	REF	QTY
Pin to Rod Coupling	MOI 55038115	8
Connecting Rod	MOI 65340250	4
Standard Bone Screws	MOI 15417120	4
Connecting Rod	MOI 15338150	2
Rod to Rod Coupling	MOI 55042011	2

Used for: Tibia Comminuted