

PROTEA

Headless Compression Screws

Introduction

The PROTEA™ Headless Compression Screw is an advanced solution designed for repairing fractures, performing joint fusions, and fixing osteotomies throughout the upper and lower extremities. Advanced implant technology and straightforward instrumentation enable the PROTEA™ systems to be an effective means of fixation for a wide variety of indications.

Implant Information

Indications


Combining versatility with an excellent rate of fixation and ongoing compression, the PROTEA™ is a great choice for a wide variety of indications, including:

- Scaphoid Fractures
- Subtalar Fusions
- Radial Head Fractures
- IP Fusions



Headless Compression Cannulated Screw - MICRO

REF (SS)	REF (TIT)	Length	Tip Diameters	Tail Diameters
MOI 16301008	MOI 36301008	8mm	2.5mm	2.8mm
MOI 16301009	MOI 36301009	9mm	2.5mm	2.8mm
MOI 16301010	MOI 36301010	10mm	2.5mm	2.8mm
MOI 16301011	MOI 36301011	11mm	2.5mm	2.8mm
MOI 16301012	MOI 36301012	12mm	2.5mm	2.8mm
MOI 16301013	MOI 36301013	13mm	2.5mm	2.8mm
MOI 16301014	MOI 36301014	14mm	2.5mm	2.8mm



Headless Compression Cannulated Screw - MINI

REF (SS)	REF (TIT)	Length	Tip Diameters	Tail Diameters
MOI 16302016	MOI 36302016	16mm	3.5mm	3.6mm
MOI 16302018	MOI 36302018	18mm	3.5mm	3.6mm
MOI 16302020	MOI 36302020	20mm	3.5mm	3.6mm
MOI 16302022	MOI 36302022	22mm	3.5mm	3.6mm
MOI 16302024	MOI 36302024	24mm	3.5mm	3.6mm
MOI 16302026	MOI 36302026	26mm	3.5mm	3.6mm
MOI 16302028	MOI 36302028	28mm	3.5mm	3.6mm
MOI 16302030	MOI 36302030	30mm	3.5mm	3.6mm



Headless Compression Cannulated Screw - STANDARD

REF (SS)	REF (TIT)	Length	Tip Diameters	Tail Diameters
MOI 16303016	MOI 36303016	16mm	4.0mm	4.1mm
MOI 16303018	MOI 36303018	18mm	4.0mm	4.1mm
MOI 16303020	MOI 36303020	20mm	4.0mm	4.1mm
MOI 16303022	MOI 36303022	22mm	4.0mm	4.1mm
MOI 16303024	MOI 36303024	24mm	4.0mm	4.1mm
MOI 16303026	MOI 36303026	26mm	4.0mm	4.1mm
MOI 16303028	MOI 36303028	28mm	4.0mm	4.1mm
MOI 16303030	MOI 36303030	30mm	4.0mm	4.1mm



USA OFFICE

4000 Island Blvd, Suite 2704
Aventura, Florida 33160

PUERTO RICO OFFICE

1660 Santa Ana Street
Esq. Fidalgo Diaz,
San Juan, Puerto Rico 00909

E-mail:

info@madisonorthoinc.com

Website:

www.madisonorthoinc.com

