


# ROCKY™ Tibia Nailing System

ROCKY™ Tibia Nailing System is suitable for all fractures extending from 7-8 cm distal to the tibial plateau, to within 5.5 cm of the distal articular surface, provided that the epiphyses are closed. ROCKY™ is a set of intramedullary nails which offers several advantages over existing systems. A major advantage of the system is the ability to insert both proximal and distal locking screws accurately and quickly without the use of X-rays, using an external mechanical targeting device.



## Features & Benefits

- Minimal invasive technique and biomechanical stability allow early functional exercise.
- Accurate and fast mechanical aiming device for proximal and distal locking.
- A locking screw of a smooth shank configuration is much stronger for a given diameter than a fully threaded screw.
- All of the implants are available in Titanium Alloy and Stainless Steel.
- The nails and locking screws are color coded for easy identification.



## Indications for Use

ROCKY™ Tibia Nailing System is indicated for: open or closed tibia shaft fractures with a very proximal and very distal extent in which locking screw fixation can be obtained, multi-fragment fractures, segmental fractures, proximal or distal non-unions, proximal or distal mal-unions, pseudarthrosis, corrective osteotomies, pathologic and impending pathologic fractures, tumor resections and comminuted fractures with or without bone loss.

## Implant Information



### End Cap, for ROCKY™ Tibia Nail

REF (S. S)	REF (TIT)	SIZE mm
MOI 11701100	MOI 31701100	0
MOI 11701000	MOI 31701000	2
MOI 11701005	MOI 31701005	5
MOI 11701010	MOI 31701010	10
MOI 11701015	MOI 31701015	15



### Locking Screw Ø4.0mm

REF (S. S)	REF (TIT)	DIA mm	LENGTH mm
MOI 11501020	MOI 31501020	4.0	20
MOI 11501022	MOI 31501022	4.0	22
MOI 11501024	MOI 31501024	4.0	24
MOI 11501026	MOI 31501026	4.0	26
MOI 11501028	MOI 31501028	4.0	28
MOI 11501030	MOI 31501030	4.0	30
MOI 11501032	MOI 31501032	4.0	32
MOI 11501034	MOI 31501034	4.0	34
MOI 11501035	MOI 31501035	4.0	35
MOI 11501036	MOI 31501036	4.0	36
MOI 11501038	MOI 31501038	4.0	38
MOI 11501040	MOI 31501040	4.0	40
MOI 11501042	MOI 31501042	4.0	42
MOI 11501044	MOI 31501044	4.0	44
MOI 11501045	MOI 31501045	4.0	45
MOI 11501046	MOI 31501046	4.0	46
MOI 11501048	MOI 31501048	4.0	48
MOI 11501050	MOI 31501050	4.0	50
MOI 11501052	MOI 31501052	4.0	52
MOI 11501054	MOI 31501054	4.0	54
MOI 11501055	MOI 31501055	4.0	55
MOI 11501056	MOI 31501056	4.0	56
MOI 11501058	MOI 31501058	4.0	58
MOI 11501060	MOI 31501060	4.0	60
MOI 11501064	MOI 31501064	4.0	64
MOI 11501065	MOI 31501065	4.0	65
MOI 11501068	MOI 31501068	4.0	68
MOI 11501070	MOI 31501070	4.0	70
MOI 11501072	MOI 31501072	4.0	72
MOI 11501076	MOI 31501076	4.0	76
MOI 11501080	MOI 31501080	4.0	80
MOI 11501085	MOI 31501085	4.0	85
MOI 11501090	MOI 31501090	4.0	90
MOI 11501100	MOI 31501100	4.0	100

### ROCKY™ Tibia Nail




REF (S. S)	REF (TIT)	DIA mm	LENGTH mm	TYPE
MOI 11301240	MOI 31301240	8.0	240	Cannulated
MOI 11301260	MOI 31301260	8.0	260	Cannulated
MOI 11301280	MOI 31301280	8.0	280	Cannulated
MOI 11301300	MOI 31301300	8.0	300	Cannulated
MOI 11301320	MOI 31301320	8.0	320	Cannulated
MOI 11301340	MOI 31301340	8.0	340	Cannulated
MOI 11301360	MOI 31301360	8.0	360	Cannulated
MOI 11301380	MOI 31301380	8.0	380	Cannulated
MOI 11301400	MOI 31301400	8.0	400	Cannulated
MOI 11302240	MOI 31302240	9.0	240	Cannulated
MOI 11302260	MOI 31302260	9.0	260	Cannulated
MOI 11302280	MOI 31302280	9.0	280	Cannulated
MOI 11302300	MOI 31302300	9.0	300	Cannulated
MOI 11302320	MOI 31302320	9.0	320	Cannulated
MOI 11302340	MOI 31302340	9.0	340	Cannulated
MOI 11302360	MOI 31302360	9.0	360	Cannulated
MOI 11302380	MOI 31302380	9.0	380	Cannulated
MOI 11302400	MOI 31302400	9.0	400	Cannulated
MOI 11303240	MOI 31303240	10.0	240	Cannulated
MOI 11303260	MOI 31303260	10.0	260	Cannulated
MOI 11303280	MOI 31303280	10.0	280	Cannulated
MOI 11303300	MOI 31303300	10.0	300	Cannulated
MOI 11303320	MOI 31303320	10.0	320	Cannulated
MOI 11303340	MOI 31303340	10.0	340	Cannulated
MOI 11303360	MOI 31303360	10.0	360	Cannulated
MOI 11303380	MOI 31303380	10.0	380	Cannulated
MOI 11303400	MOI 31303400	10.0	400	Cannulated
MOI 11391240	MOI 31391240	11.0	240	Cannulated
MOI 11391260	MOI 31391260	11.0	260	Cannulated
MOI 11391280	MOI 31391280	11.0	280	Cannulated
MOI 11391300	MOI 31391300	11.0	300	Cannulated
MOI 11391320	MOI 31391320	11.0	320	Cannulated
MOI 11391340	MOI 31391340	11.0	340	Cannulated
MOI 11391360	MOI 31391360	11.0	360	Cannulated
MOI 11391380	MOI 31391380	11.0	380	Cannulated
MOI 11391400	MOI 31391400	11.0	400	Cannulated
MOI 11392240	MOI 31392240	12.0	240	Cannulated
MOI 11392260	MOI 31392260	12.0	260	Cannulated
MOI 11392280	MOI 31392280	12.0	280	Cannulated
MOI 11392300	MOI 31392300	12.0	300	Cannulated
MOI 11392320	MOI 31392320	12.0	320	Cannulated
MOI 11392340	MOI 31392340	12.0	340	Cannulated
MOI 11392360	MOI 31392360	12.0	360	Cannulated
MOI 11392380	MOI 31392380	12.0	380	Cannulated
MOI 11392400	MOI 31392400	12.0	400	Cannulated

### ROCKY™ Tibia Nail


REF (S. S)	REF (TIT)	DIA mm	LENGTH mm	TYPE
MOI 11001240	MOI 31001240	8.0	240	Solid
MOI 11001260	MOI 31001260	8.0	260	Solid
MOI 11001280	MOI 31001280	8.0	280	Solid
MOI 11001300	MOI 31001300	8.0	300	Solid
MOI 11001320	MOI 31001320	8.0	320	Solid
MOI 11001340	MOI 31001340	8.0	340	Solid
MOI 11001360	MOI 31001360	8.0	360	Solid
MOI 11001380	MOI 31001380	8.0	380	Solid
MOI 11001400	MOI 31001400	8.0	400	Solid
MOI 11002240	MOI 31002240	9.0	240	Solid
MOI 11002260	MOI 31002260	9.0	260	Solid
MOI 11002280	MOI 31002280	9.0	280	Solid
MOI 11002300	MOI 31002300	9.0	300	Solid
MOI 11002320	MOI 31002320	9.0	320	Solid
MOI 11002340	MOI 31002340	9.0	340	Solid
MOI 11002360	MOI 31002360	9.0	360	Solid
MOI 11002380	MOI 31002380	9.0	380	Solid
MOI 11002400	MOI 31002400	9.0	400	Solid
MOI 11003240	MOI 31003240	10.0	240	Solid
MOI 11003260	MOI 31003260	10.0	260	Solid
MOI 11003280	MOI 31003280	10.0	280	Solid
MOI 11003300	MOI 31003300	10.0	300	Solid
MOI 11003320	MOI 31003320	10.0	320	Solid
MOI 11003340	MOI 31003340	10.0	340	Solid
MOI 11003360	MOI 31003360	10.0	360	Solid
MOI 11003380	MOI 31003380	10.0	380	Solid
MOI 11003400	MOI 31003400	10.0	400	Solid
MOI 11070240	MOI 31070240	11.0	240	Solid
MOI 11070260	MOI 31070260	11.0	260	Solid
MOI 11070280	MOI 31070280	11.0	280	Solid
MOI 11070300	MOI 31070300	11.0	300	Solid
MOI 11070320	MOI 31070320	11.0	320	Solid
MOI 11070340	MOI 31070340	11.0	340	Solid
MOI 11070360	MOI 31070360	11.0	360	Solid
MOI 11070380	MOI 31070380	11.0	380	Solid
MOI 11070400	MOI 31070400	11.0	400	Solid
MOI 11071240	MOI 31071240	12.0	240	Solid
MOI 11071260	MOI 31071260	12.0	260	Solid
MOI 11071280	MOI 31071280	12.0	280	Solid
MOI 11071300	MOI 31071300	12.0	300	Solid
MOI 11071320	MOI 31071320	12.0	320	Solid
MOI 11071340	MOI 31071340	12.0	340	Solid
MOI 11071360	MOI 31071360	12.0	360	Solid
MOI 11071380	MOI 31071380	12.0	380	Solid
MOI 11071400	MOI 31071400	12.0	400	Solid

## INSTRUMENTS INFORMATION

### ROCKY Tibia Nail Instrument



REF	DESCRIPTION	UNITS
1	MOI 0201014 Rigid Reamer, 11.0mm	1
2	MOI 0201015 Rigid Reamer, 10.0mm	1
3	MOI 0201016 Rigid Reamer, 9.0mm	1
4	MOI 0100003 Drill Bit 4.0mm, with Drill Stop	2
5	MOI 0506001 T-Handled Reamer 4.0mm, for Position Rod	1
6	MOI 0201007 T-Handled Screwdriver SW3.5	1
7	MOI 0201008 Depth Gauge for Locking Screw	1
8	MOI 0201003 Point Awl	1
9	MOI 0100023 Allen Wrench SW5, Large	1
10	MOI 0201005 Soft Tissue Protector, 8.0mm	1
11	MOI 0100025 Allen Wrench SW3, Small	1
12	MOI 0506002 Protection Sleeve 10.0/8.0mm, for Position Rod	1
13	MOI 0201010 Quick Coupling T-Handle	1
14	MOI 0201002 Protection Sleeve 8.0/4.1mm, for Locking Screw	1
15	MOI 0201001 Hex Wrench, Open	1
28	MOI 0100043 Flexible Reamer Shaft	1
32	MOI 0200003 Replaceable Reamer Heads	3



REF	DESCRIPTION	UNITS
16	MOI 0506003 Insertion Handle	1
17	MOI 0506004 Distal Aiming Arm	1
18	MOI 0506005 Proximal Aiming Arm	1
19	MOI 0201031 Handle, for Nail Extraction	1
20	MOI 0201029 L-Position Rod	1
21	MOI 0506006 Connecting Bolt	1
22	MOI 0201025 Locking Screw Extractor	1
23	MOI 0506007 U-Shape Block 8/9/10, for Stabilizing Rod	3
24	MOI 0506008 Position Rod	1
25	MOI 0201023 T-Handled Tap 6.0mm	1
26	MOI 0100040 Sliding Hammer, for Nail Insertion / Extraction	1
27	MOI 0506009 Guide Bar	1
29	MOI 0100044 Guide Wire, with Ball Tip	1
30	MOI WR00001 Locking Wheel, Solid, for Rocky Guide Bar	1
31	MOI WR00002 Locking Wheel, Cannulated	1

REF	DESCRIPTION	UNITS
33	MOI 0506100 Tibia Nail Container	



### **USA OFFICE**

4000 Island Blvd, Suite 2704  
Aventura, Florida 33160

### **PUERTO RICO OFFICE**

1660 Santa Ana Street  
Esq. Fidalgo Diaz,  
San Juan, Puerto Rico 00909

E-mail:

[info@madisonorthoinc.com](mailto:info@madisonorthoinc.com)

Website:

[www.madisonorthoinc.com](http://www.madisonorthoinc.com)

